




User Guide for Lingolet Interpreter Video

Connect to an interpreter with the Lingolet Interpreter Video in a few simple steps:

1. Power on the device and connect to wi-fi:



2. Select the language that you need to be interpreted on the home screen:

Which language do you need an interpreter for?

 Spanish	 Chinese	 Japanese	 Russian
 Ukrainian	 American Sign Language	 Vietnamese	 Settings


Interpreter Gender: No preference

Industry: Medical


 Audio Call

 Video Call


3. Click either **Audio Call** or **Video Call** to connect to an interpreter:



4. When you are finished with the call, rate the interpreter and call quality:



5. Click **Settings** on the home screen to change languages, your industry, or the gender of interpreters:



Add New Languages to the Home Screen

1. Click on the **Settings** button.
2. Click on the **Target Language** dropdown and select a new language.
3. The last language in the list under **Target Languages on Home Screen** will be removed, and the new language will be added to the list.

Changing Industry or Gender Preference

1. Click on the **Settings** button.
2. Click on the **Gender** dropdown and select the interpreter's gender.
3. Click on the **Industry** dropdown and select your industry.

If you have any questions, contact Lingolet Support at support@lingolet.com or [Lingolet.com/support/](https://lingolet.com/support/)

Copyright Notice

This document is © 2023 by Lingolet, Inc. USA. This document is protected by copyright. All rights reserved. No part of this document may be reproduced in any form by any means without prior written permission from Lingolet, Inc.